

Enhanced HiRise Effective Doses⁽¹⁾

	Lite Protocol mSv	Std Protocol mSv	Lrg Protocol mSv
HIP	1.9	3.1	2.5
kV	120	120	120
mA	20	20	20
Exposure Time (Sec)	9.36	14.4	25.2
PELVIS	Not Calculated	4.3	4.5
kV	120	120	120
mA	20	20	20
Exposure Time (Sec)	9.36	14.4	25.2
KNEE*	0.0117	0.0178	0.0214
kV	100	120	120
mA	12	14	18
Exposure Time (Sec)	10.8	10.8	10.8

	Lite Protocol mSv	Std Protocol mSv	Lrg Protocol mSv
HAND/WRIST	0.0024	0.0056	0.0086
kV	100	120	120
mA	12	12	12
Exposure Time (Sec)	5.76	5.76	6.24
FOOT/ANKLE (LFOV)*	0.0108	0.0181	0.0293
kV	100	120	120
mA	12	12	20
Exposure Time (Sec)	8.64	8.64	9.36
FOOT/ANKLE (MFOV)*	0.0132	0.0215	0.0424
kV	100	120	120
mA	12	12	20
Exposure Time (Sec)	8.64	8.64	9.36

	Lite Protocol mSv	Std Protocol mSv	Lrg Protocol mSv
PEDIATRIC KNEE*			
5 Years Old	0.0361	0.0262	0.0139
10 Years Old	0.0262	0.0430	0.0217
15 Years Old	0.0721	0.0511	0.0260
PEDIATRIC HAND/WRIST			
5 Years Old	0.0075	0.0193	0.0307
10 Years Old	0.0043	0.0108	0.0169
15 Years Old	0.0024	0.0057	0.0086
PEDIATRIC FOOT/ANKLE (LFOV)*			
5 Years Old	0.0314	0.0632	0.1030
10 Years Old	0.0163	0.0294	0.0478
15 Years Old	0.0104	0.0174	0.0282
PEDIATRIC FOOT/ANKLE (MFOV)*			
5 Years Old	0.0394	0.0751	0.1460
10 Years Old	0.0200	0.0346	0.0681
15 Years Old	0.0127	0.0206	0.0407

VARIOUS DIAGNOSTIC RADIOLOGY PROCEDURES, INCLUDING FLOUROSCOPY⁽²⁾

	mSv
EXTREMITY	<0.001
HIP	0.7
PELVIS	0.6
MEDICAL CT	
	mSv
UNILATERAL FOOT & ANKLE	0.07 ⁽³⁾
STANDARD HIP	3.5 ⁽¹⁾
LOW DOSE HIP	0.97 ⁽⁴⁾

kV, mA, and Exposure Times are the same as those for adult protocols.

Sources:

* Bilateral Exposure.
 (1) Internal study overseen by Jaydev Dave, PhD, DABR, MS, FAAP. Abstract submitted for review.
 (2) Mettler FA, et al: Effective Doses in Radiology and Diagnostic NuclearMedicine: A Catalog, Radiology 2008 248:254-263
 (3) Biswas D, Bible JE, Bohan M, Simpson AK, Whang PG, Grauer JN. Radiation exposure from musculoskeletal computerized tomographic scans. J Bone Joint Surg Am. 2009 Aug;91(8):1882-9. doi: 10.2106/JBJS.H.01199. PMID: 19651945.
 (4) Arthroscopyjournal.org (2019, May 01). Low-Dose Computed Tomography... Retrieved Oct. 20, 2022, from https://www.arthroscopyjournal.org/article/S0749-8063(18)31067-3/fulltext